

PAJUNK®

SonoTAP Dual Catheter Kit
for continuous techniques



Abdominal blocks


SonoTAP Dual Catheter Kit

Designed by Anaesthetist for Anaesthetists, The SonoTAP Dual catheter kit offers you an extra long echogenic needle and dual catheters with multiple lateral eyes for continuous TAP blocks procedures.

Tuohy Sono Cannula

The visible cannula for added safety


With the TuohySono, PAJUNK has risen to the challenge of market development and is the first company offering a cannula perfectly matched to the requirements placed on ultrasound-supported anaesthesia.



EpiLong Catheter

maximum stability and flexibility in a catheter

The new EpiLong catheter developed by PAJUNK to combines two essential attributes for atraumatic and safe puncturing. It is absolutely soft while being absolutely stable at the same time. Six openings on the distal catheter end facilitate perfectly uniform distribution of the anaesthetic.



The six lateral openings on the catheter tip guarantees perfectly uniform distribution of the anaesthetic, which provides optimal development of the effectiveness.

Ordering information

SonoTAP Dual Catheter Kit



Item No.	Product	PU
101151-603	SonoTAP Dual Catheter Kit 1 x TuohySono cannula with "Cornerstone" reflectors 18G x 150mm 2 x Catheter 20G x 90 cm, with 6 lateral openings on the first 18 mm of the catheter tip 2 x Clamping Adapters (Red) 2 x Filters	10/Box

For more information about this product please don't hesitate to contact us



www.surimex.com.au

SURIMEX

Head Office

Pajunk: Australia and New Zealand

15 Redland Drive Mitcham

VIC 3132 Australia

T +613 8833 3804

F +613 8833 3809

PAJUNK GmbH

Medizintechnologie

Karl-Hall-Strasse 1

D-78187 Geisingen

Germany

Telefon +49 (0) 77 04/92 91-0

Telefax +49 (0) 77 04/92 91-6 00

www.pajunk.com