

# SURIMEX

## ***Silicone sheeting & Silicone tubing***

*Surimex offers various sizes of sterile and non sterile silicone products.*



### ***Quality Specification***

Type	Silicone Rubber, ST-EC-60-001
Hardness	60+/-5D. Shore, DIN 53 505
Density	1.17 g/ccm, DIN 53 479A
Colour	Translucent (or as required)

### ***Characteristics***

- Heat resistant in hot air up to 200°C
- Suitable for autoclaving (135°C= 2.2-2.5 bar)
- Suitable for gamma- sterilization
- Suitable for application in the medical technology and food areas

## Silicone Sheeting



Silicone Sheeting can be used in reconstructive surgery. It is available with or without reinforcing and in various thicknesses.

## Silicone Sheeting

Silicone Sheeting uses include: nasal splinting, wound dressing, wound drains, abdominal vents or coverings for neonatal gastroschisis, external suture bolsters, laboratory uses and other uses

## Reinforced Silicone Sheeting

Reinforced silicone sheet (0.4mm) uses include: osteogenesis, repair of urethral structures and duramater, repair of orbital floor fractures (Special Order).

## Ordering Information

Silicone Sheeting: non sterile

REF Code	Thickness (mm)	Length (mm)	Width (mm)	U/M
245157	1	300	275	1
5101002	2	300	275	4*

(\*Min. Order Qty)

Reinforced Silicone Sheeting: non sterile

REF Code	Thickness (mm)	Length (mm)	Width (mm)	U/M
TP1701-15	0.4	150	200	1

## Silicone Tubing



## Silicone Tubing: sterile

Sterile silicone tubing is pre packaged (sterile)

## Silicone Tubing: non sterile

Non sterile silicone tubing is available in coils with various inner diameters (I.D.) and outer diameters (O.D.)

## Ordering Information

General Purpose Silicone Tubing: Sterile

























REF Code	I.D. (mm)	O.D. (mm)	Length (cm)	U/M
246045-P	5	8	80	40/ Box

Silicone Snugger 3/Pack: Sterile

REF Code	FG	I.D. (mm)	O.D. (mm)	Length (cm)	U/M
222253-P	14	5	7	15	50/ Box

## Ordering Information

## Silicone Tubing: non-sterile

REF Code	I.D (mm)	O.D (mm)	Wall Width	Coil Length (mt)	100% Size	REF Code	I.D (mm)	O.D (mm)	Wall Width	Coil Length (mt)	100% Size
246001	1	2	0.5	25							
246002	1	3	1	25		246026	9	12	1.5	25	
246003	2	4	1	25							
246004	3	5	1	25		246028	10	14	2	25	
246005	3	6	1.5	25							
246006	3	7	2	25		246031	10	15	2.5	25	
246007	4	6	1	25							
246008	4	7	1.5	25							
246009	4	8	2	25							
246010	5	7	1	25							
246011	5	8	1.5	25							
246011G											
Coloured Silicone Tubing (Green)											
246012	5	9	2	25							
246014	6	8	1	25							
246015	6	9	1.5	25							
246016	6	10	2	25							
246017	6	11	2.5	25							
246019	7	10	1.5	25							
246020	7	11	2	25							
246023	8	11	1.5	25							
246024	8	12	2	25							
246025	8	14	3	25							

For all stock availability please head to our website

[www.surimex.com.au](http://www.surimex.com.au) or call +613 8833 3804

*Do you have any questions? Our customer service & sales support will be glad to assist you.*

Products Distributed in Australia and New Zealand by

**SURIMEX**

15 Redland Drive  
Mitcham VIC 3132  
Australia

T (613) 8833 3804  
F (613) 8833 3809  
E [sales@surimex.com.au](mailto:sales@surimex.com.au)

[www.surimex.com.au](http://www.surimex.com.au)