

## Extension Lines

---

extension tubing for intravenous administration



# Extension tubes

Extension tubing for one-way or multi-way intravenous infusion is a sterile, single-use medical device.

Extension tubing is used for the administration of drug substances or sample collections from veins or solely for sample collections from arteries in combination with other infusion devices: infusion pumps, manifolds, stopcocks, bidirectional valves, etc.

Extension tubing for intravenous infusion helps administer the medicinal product by carrying the drug substances into the patient's bloodstream via a catheter. Extension tubing is also used to collect blood.

Compatible with lipid solutions, chemotherapy solutions, blood and blood derivatives.



Fixed Male to Female Luer Lock connections



DEHP (DOP), latex, BBP, DBP, DIBP or DHP Free

Item no.	Description	I.D (mm)	O.D (mm)	Length (cm)	Connection type	Tubing Material	PU	Stock availability
PB3120	Extension Line	1.0	2.5	200	Female / Male	PE/PVC	40	stocked item
PB3202	Extension Line	1.5	3.0	25	Female / Male	PE/PVC	40	stocked item
PN3220	Extension Line	1.5	3.0	200	Female / Male	PVC	40	stocked item
PB1220	Extension Line	1.5	3.0	200	Male / Male	PE/PVC	40	stocked item

MD class: IIa CE0459 (Please read the instruction for use.)

*Do you have any questions? Our customer service & sales support will be glad to assist you.*

Products distributed in Australia and New Zealand by

**SURIMEX**

15 Redland Drive  
Mitcham VIC 3132  
Australia

(613) 8833 3804  
sales@surimex.com.au  
[www.surimex.com.au](http://www.surimex.com.au)