


SURIMEX

Comvac mini

Compression vacuum system



Comvac mini compression vacuum system

A manual compression system for post operative wound drainage. They provide a constant low negative pressure during wound healing process to prevent the wound surfaces and underlying layer from being disturbed.

Comvac mini bellow would typically be used after: Minor surgeries, such as head, neck, hand, face and plastic surgery.



Compression bellow

- Easy to compress bellow




Screw connection

- Screw connection for use of use
- User friendly




Connection holes

- Offering connection points for users if required



Drain with trocar

- Pre-connected trocar for user



Drain depth markers

- Depth markers at 1cm & 2cm from beginning of drain perforations

Item no.	Description	Volume	Vacuum	Drain	PU	Stock availability
220210	Comvac mini set	50ml	Low	6	10	stocked item

Do you have any questions? Our customer service & sales support will be glad to assist you.

Products distributed in Australia and New Zealand by

SURIMEX

15 Redland Drive
Mitcham VIC 3132
Australia

T (613) 8833 3804
E sales@surimex.com.au
www.surimex.com.au