


SURIMEX

Manifolds, Universal spike adapter, Scissor clamps




Manifolds, Universal spike adapter, Scissor clamps

Manifolds are made of poly carbonate. Each manifold is a one-piece moulding production and has a holding plate for easy and secure attachment to the IV pole with a special clamp. They are transparent and pressure resistant till 2 bar. The stopcock handles are blue or with different colours to differentiate the monitoring lines.

Universal Spike Adapter The exchange coupler can be used to transfer liquids from vials into syringes etc.

Scissor Clamp Plastic scissor clamps can be used for clamping tubes without damage.


Manifolds



Universal spike adapter



Scissor clamps



Item no.	Description	PU	Stock availability
701030	Manifold (3 port) (luer lock, 1 male & 4 female)	50	stocked item
701050	Manifold (5 port) (luer lock, 1 male & 6 female)	50	stocked item
210817-P	Universal Spike Adapter: Length 75mm (spike / Luer lock end)	100	stocked item
247001-P	Scissor clamps (blue) (non-sterile)	50	stocked item

Do you have any questions? Our customer service & sales support will be glad to assist you.

Products distributed in Australia and New Zealand by

SURIMEX

15 Redland Drive
Mitcham VIC 3132
Australia

T (613) 8833 3804
E sales@surimex.com.au
www.surimex.com.au