


SURIMEX

Safety Huber Needle

Power injection, non coring, safety huber needle




Safety Huber Needle

Elevate your medical procedures with our Safety Huber Needle. Crafted to endure pressures up to 325psi, this needle excels in power injection of contrast media. Its non-coring design ensures minimal tissue disruption, prioritizing patient comfort. Experience the future of safety and precision in medical practices – experience the Safety Huber Needle.



Needle safety design

- Audible click
- Preventing needle stick injuries



Huber needle

- Non coring needle
- Minimising damage to port septum




Colour coded wings

- Needle gauge identification



Y-site with needlessly connection

- Offers a pre connected safety
- Transparent for easy visualisation
- Needleless self closing septum
- Easy cleaning and disinfection



Transparent base

- Optimal site visualisation



Wide base with padding

- Needle stability
- Extra padding for patient comfort

Item no.	Description	Colour	PSI Rating	PU	Stock availability
VA032020SP	Safety Huber Needle: 20G x 19mm	Yellow	325	50	stocked item
VA032025S	Safety Huber Needle: 20G x 25mm	Yellow	325	50	available in future
VA032038S	Safety Huber Needle: 20G x 38mm	Yellow	325	50	available in future
VA032019SY	Safety Huber Needle: 20G x 19mm with needlessly y-site	Yellow	325	50	available in future
VA032025SY	Safety Huber Needle: 20G x 25mm with needlessly y-site	Yellow	325	50	available in future

* Please note other sizes available upon request

Do you have any questions? Our customer service & sales support will be glad to assist you.

Products distributed in Australia and New Zealand by

SURIMEX

15 Redland Drive
Mitcham VIC 3132
Australia

T (613) 8833 3804
E sales@surimex.com.au
www.surimex.com.au