


PAJUNK®



InfiltraLong System

Continuous Wound Infiltration



InfiltraLong (with split cannula)



Item no.	needle	InfiltraLong Catheter	Kit includes	PU
001158-00Z	split cannula with puncture needle	19G x 420mm 15 openings along the first 40mm with helical coil	Filter FixoLong	10
001158-20B	split cannula with puncture needle	19G x 600mm 60 openings along the first 150mm with helical coil	Filter FixoLong	10
001158-20D	split cannula with puncture needle	19G x 600mm 45 openings along the first 100mm with helical coil	Filter FixoLong	10
001158-30C	split cannula with puncture needle	19G x 700mm 88 openings along the first 220mm with helical coil	Filter FixoLong	10

InfiltraLong (with Tuohy needle)



Item no.	Needle	InfiltraLong Catheter	Kit includes	PU
001159-20D	Tuohy 17G x 90mm	19G x 600mm 15 openings along the first 40mm with helical coil	Clamp adaptor Filter & FixoLong	10

InfiltraLong Sono (Tuohy Sono needle)



Item no.	Needle	InfiltraLong Catheter	Kit includes	PU
001185-20B	Tuohy Sono 17G x 90mm	19G x 600mm 15 openings along the first 40mm with helical coil	Clamp adaptor Filter & FixoLong	10

Y-connector tube



Item no.	Item no.	PU
001151-37W	Y-connector tube Injecting tube with y-connection for the simultaneous use of 2 catheters (Luer lock ends)	10

FixoCath



Item no.	Item no.	PU
001161-37Z	FixoCath Catheter fixation device for patient comfort	10

PAJUNK® GmbH Medizintechnologie
Global Headquarters
Karl-Hall-Strasse 1
78187 Geisingen, Germany
Tel.: +49 7704 9291-0
Fax: +49 7704 9291-600
info@pajunk.com
pajunk.com

Do you have any questions? Our customer service & sales support will be glad to assist you.

Products distributed in Australia and New Zealand by

SURIMEX

15 Redland Drive
Mitcham VIC 3132
Australia

T (613) 8833 3804
E sales@surimex.com.au

www.surimex.com.au