


# SURIMEX

## Safety Huber Needle

---

Power injection, non coring, safety huber needle




# Safety Huber Needle

Elevate your medical procedures with our Safety Huber Needle. Crafted to endure pressures up to 325psi, this needle excels in power injection of contrast media. Its non-coring design ensures minimal tissue disruption, prioritizing patient comfort. Experience the future of safety and precision in medical practices – experience the Safety Huber Needle.



Needle safety design

- Audible click
- Preventing needle stick injuries



Huber needle

- Non coring needle
- Minimising damage to port septum




Colour coded wings

- Needle gauge identification



Printed Clamp: PSI pressure

- Quick and easy Identification
- PSI rating & max. mls per second
  - CT Identification stickers: Still available upon request



Transparent base

- Optimal site visualisation



Wide base with padding

- Needle stability
- Extra padding for patient comfort

Item no.	Description	Colour	PSI Rating	PU	Stock availability
VA032015S	Safety Huber Needle: 20G x 15mm	Yellow	325	50	stocked item
VA032020SP	Safety Huber Needle: 20G x 19mm	Yellow	325	50	stocked item
VA032025S	Safety Huber Needle: 20G x 25mm	Yellow	325	50	stocked item
VA032038S	Safety Huber Needle: 20G x 38mm	Yellow	325	50	available in future
VA032019SY	Safety Huber Needle: 20G x 19mm with needleless y-site	Yellow	325	50	available in future
VA032025SY	Safety Huber Needle: 20G x 25mm with needleless y-site	Yellow	325	50	available in future

\* Please note other sizes available upon request

*Do you have any questions? Our customer service & sales support will be glad to assist you.*

Products distributed in Australia and New Zealand by



15 Redland Drive  
Mitcham VIC 3132  
Australia

T (613) 8833 3804  
E sales@surimex.com.au  
[www.surimex.com.au](http://www.surimex.com.au)